

Compound	Molecular Mass (g/ mol)	Liquid AN Medium		
		unmodified	N-limited	P-limited
Nanopure water			1000 mL	
Base salts ^a				
NaCl	58.4		370 mM	
MgSO ₄ · 7H ₂ O	246.5		30 mM	
MgCl ₂ · 6H ₂ O	203.3		20 mM	
CaCl ₂ · 2H ₂ O	147.0		10 mM	
KCl	74.6		10 mM	
Nutrients ^b				
EDTA (disodium salt)	292.2		17 µM	
Na ₂ CO ₃ · H ₂ O	124.0		161 µM	
NaHCO ₃	84.0		500 µM	
NaNO ₃	85.0	176 µM	220 µM	440 µM
K ₂ HPO ₄	174.2	88 µM	22 µM	11 µM
Cyano trace metals ^c			1 mL	
ZnSO ₄ · 7H ₂ O (0.222 g/L)				
MnCl ₂ · 4H ₂ O (1.4 g/L)				
Co(NO ₃) ₂ · 6H ₂ O (0.025 g/L)				
Na ₂ MoO ₄ · 2H ₂ O (0.39 g/L)				
Citric Acid · H ₂ O (6.25 g/L)				
Ferric Ammonium Citrate (6 g/L)				
V _A Vitamin Solution ^d			1 mL	
	<i>Stock Concentration</i>	<i>To prepare solution</i>		
Inositol		100 mg		
Thiamine · HCl		20 mg		
Vitamin B ₁₂	1 g/L	0.1 mL		
Biotin	0.1 g/L	1 mL		
Folic Acid	2 g/L	0.1 mL		
p-aminobenzoic acid	2 g/L	0.5 mL		
Niacin (Nicotinic acid)	1 g/L	10 mL		
Ca d-pantothenate	2 g/L	10 mL		
Pyridoxine	1 g/L	10 mL		

^a Combine all salts, bring to final solution volume of 1 liter then autoclave

^b Add sterile stock solution to autoclaved salts

^c To prepare stock, dissolve each separately and bring to 1 liter

^d V_A vitamin solution is added for enrichment and purification; to prepare stock, prepare separate stock solutions, then add to 60 mL nanopure water. Filter-sterilize through 0.1 µm filter and store at -20°C in 50 mL aliquots.